



Ultrasonic Tooth Cleaner Fairywill FW-C8 (black)

Ref: 6973734202658

Ultrasonic Tooth Cleaner Fairywill FW-C8 (black)

FairyWill FW-C8 ultrasonic tartar remover (black)

Get rid of even extremely stubborn tartar with the FairyWill ultrasonic device. The FW-C8 model supports 3 modes of operation to meet different user needs. It features IPX7 water resistance, is equipped with a smart switch, and one charge is enough to enjoy its operation for up to 30 days. It comes with interchangeable tips and a mirror to help you see into hard-to-reach areas of your mouth, so you can carry out the cleaning process more effectively. Operating the product is simple, and a light indicator will let you know the battery charge status.

Effective performance

You no longer need to use professional services to ensure proper oral hygiene, get a beautiful smile and eliminate tartar. The FW-C8 will work perfectly as your home dentist. With the help of the included tips, you will say goodbye to stains and plaque for good and regain your confidence!

3 work modes

You can easily adjust the operation of the FairyWill ultrasonic device to your needs. FW-C8 supports up to 3 modes of operation. Soft will be a great choice for people who have particularly sensitive teeth. Normal mode will help maintain proper oral hygiene. Strong, on the

other hand, will work well for the most demanding users who care about deep cleaning.

Convenience of use

Operating the FairyWill appliance will not give you any trouble. You will easily install the tip of your choice. Then, all you have to do is hold down the corresponding button for 3 seconds to start the FW-C8 and select the mode that suits you. For added convenience while cleaning your teeth, you can use the included mirror. The recommended usage time of the product is 20 to 30 minutes per day.

Long working time

One charge of the FairyWill device is enough to use it for up to 30 days. It takes about 4 hours to renew the FW-C8's energy with a USB cable. You can check the battery status with the light indicators. Red indicates that the device is close to being discharged or has been plugged in for charging. Yellow, on the other hand, is equivalent to the fact that the battery energy has been completely renewed.

Thoughtful solutions

To enhance safety, the FW-C8 automatically shuts off 10 minutes after starting operation. The solution used helps minimize the risk of overheating. In addition, the product provides IPX7 water resistance and features a non-slip design.

Included:

Ultrasonic tartar remover x1

Flat tip x2

Sharp tip x2

Mirror x1

Instruction manual x1

Manufacturer

FairyWill

Model

FW-C8

Working temperature

from -5°C to 40C

Humidity

≤80%

Operating voltage

AC12.0V

Input voltage

AC100V-240V

Rated power

0,5-3W

Range of tip movement

< 100um

Frequency of vibration

30-40kHz

Degree of water resistance

IPX7

Color

Black

Price:

€ 27.50

Health & Beauty, Sonic toothbrush