



## Cavitation Peeling with Light Therapy ANLAN ALCPJ05-02 (Whit

Cavitation Peeling with Light Therapy ANLAN ALCPJ05-02 (Whit

Ref: 6953156299184

Cavitation peeling with light therapy ANLAN ALCPJ05-02

Anlan device will allow you to take care of your skin in a quick and easy way. The peeling uses red and blue light therapy and positive and negative ions. It reduces the appearance of wrinkles, treats acne and firms the skin. It also makes it easier to get rid of blackheads and whiteheads. What's more, its operation is simple, so it won't give you any trouble.

Thoroughly cleanse the face skin

With the help of Anlan cavitation peeling you will thoroughly cleanse the skin of your face. The device uses ultrasonic vibrations and works with the amount of 24000 vibrations per second. What is more, the positive ions eliminate impurities located in the deep layers of the skin. As a result, your skin is clean and ready for your next skin care treatment.

Effective reduction of wrinkles

Restore your skin's smoothness and elasticity. Microcurrent EMS guarantees a lifting effect thanks to stimulation of facial muscles. What is more, it activates collagen production and improves skin condition. As a result the complexion is firm and wrinkles are smoothed out.

The power of negative ions

Cavitation peeling from Anlan will increase the effectiveness of the cosmetics you use. Negative ions together with high frequency

vibrations help nutrients reach deeper layers of the skin. Now you can enjoy perfectly moisturized skin.

#### Red and blue light therapy

The device also offers red and blue light therapy. Red light is responsible for increased skin elasticity and slowing down the aging process. Blue light, on the other hand, cleanses the skin of bacteria and supports acne treatment.

#### Easy to use

Operating the device is simple and convenient. Wash your face and then turn on Ultrasonic mode. After about 2-3 minutes, select the Ion+ mode for extra facial cleansing and exfoliation of dead skin. Then apply your favorite cream and activate Ion mode. Finally, use the EMS mode, for smoothing the skin.

#### Included

cavitation peeling

charging cable

user manual

#### Brand

ANLAN

#### Model

JD-CP001

#### Red light wavelength

605-700 nm

#### Blue light wavelength

450-480 nm

#### Dimensions

approx. 170 x 43 x 19 mm

#### Weight

approx. 88 g

#### Material

ABS + stainless steel

#### Rated power

approx. 1,5 W

#### Vibration frequency

24,000 vibrations per second

#### Battery capacity

500 mAh

#### Charging time

approx. 2 h

#### Automatic switch off

after approx. 10 minutes

Preço:

**€ 27.00**

Saúde e beleza, Skincare equipment