



Ultrasonic facial massager with light therapy ANLAN 01-ADRY1

Ref: 6953156301269

Ultrasonic facial massager with light therapy ANLAN 01-ADRY1

Ultrasonic facial massager with light therapy ANLAN 01-ADRY15-001

Do you want to take care of your skin, but still lack the time for detailed care? With the ultrasonic face massager from ANLAN, five minutes a day is enough to properly take care of your complexion and get satisfactory results. Operation of the device is incredibly simple and you can easily adjust the work of the massager to your needs, and you can choose from up to 4 care modes, 3 types of phototherapy and 3 levels of temperature. Thanks to the powerful battery, the device will serve you for a long time.

Effective facial cleansing

If you are struggling with acne, excess visible pores or your skin has a grayish tint, think about changing your cleansing method. When washing your face by hand doesn't bring the desired results, take advantage of the possibilities offered by ANLAN's ultrasonic facial massager. A warm compress (temperature of 38°C / 40°C / 42°C) will effectively open the pores and improve the condition of your skin. On the other hand, a cold compress (temperature of 6°C / 8°C / 10°C) will narrow the pores, make the skin more elastic and increase the effectiveness of the skin care cosmetics used.

Light therapy

Feel like you're in a spa and perform phototherapy yourself. ANLAN ultrasonic massager provides as many as 3 types of light therapy. If you want to make your skin more flexible and eliminate wrinkles, use the red light. For effective acne control, maximum skin cleansing



and elimination of bacteria, use the blue light. The last of the lights - green - improves microcirculation, eliminates puffiness and reduces sebum secretion.

Skin lifting at home

Are you noticing wrinkles appearing more and more often? Or perhaps your skin is weakened and less elastic? You don't have to make an appointment for expensive treatments to solve these problems. The EMS micro-currents that the ANLAN massager uses penetrate deep into your muscles to help make your skin more elastic and regenerate. After just 7 days of regular use (2x a day, 10 to 15 minutes) you will enjoy a healthy and youthful appearance!

One device, many uses

ANLAN's ultrasonic facial massager has so many functions that it replaces as many as 5 different skin care devices. You can use the product to perform phototherapy, deep cleanse your skin, make compresses or even slim your face. And these are still not all the capabilities of the device! The mode of operation using ions and high-frequency vibrations will ensure that the nutrients of the skin care products used will reach the deepest layers of the skin, so you will enjoy a flawless complexion. Remember to wash off your makeup thoroughly before using the massager.

Thoughtful design

The LED display will inform you about the working status of the ANLAN device. You will easily check the temperature level or the battery charge, among other things. The time needed to charge the massager is about 3 hours. Using the device will not cause you any trouble. In addition, the massager fits perfectly in your hand, so you will not feel any discomfort while taking care of your skin. A special stand is included, which allows you to store the product conveniently and safely. Every effort has been made to make the use of the ANLAN massager effective and comfortable.

Included:

Ultrasonic facial massager with light therapy

Cosmetic pad ring

Base

Cosmetic pad

USB-C charging cable

User manual

Manufacturer

ANLAN

Model

01-ADRY15-001

Material

ABS + stainless steel

Dimensions

179 x 60 x 61 mm

Weight

195 g

Vibration frequency

8000 RPM

Light wavelength

Red light 622 nm Blue light 467 nm Green light 523 nm

Battery capacity

1000 mAh

Charging time

3 h



Usage time

approximately 62 minutes

Automatically power off

Clean, EMS, Moist mode: 5 minutes Cool mode: 2 minutes

Preço:

Antes: € 55.0056

Agora: € 54.50

Saúde e beleza, Health and beauty, Skincare equipment, Facial care devices, Massagers

