



Body slimming device ANLAN 02-AXTY31-02A
Body slimming device ANLAN 02-AXTY31-02A

Slimming body device ANLAN 02-AXTY31-02A

Want to improve the appearance of the skin on your thighs, abdomen and arms? Regular manual care will help you achieve the desired effect without much sacrifice. This massage device offers 3 modes to comprehensively care for your body. Indulge yourself with a pleasure that will have the side effect of improving your metabolism and smoothing and improving the elasticity of your skin!

Ref: 6953156301559

Telefone: +351962484153

Email: info@mobilemax.pt

(Chamada para rede móvel nacional)

Blue light power

Mode one uses ultrasonic waves in combination with vibration and blue light. The vibration motor generates up to 9,000 rpm, providing a pleasant massage experience. This configuration inhibits sebum secretion, acting as an anti-inflammatory and antiseptic. Mode one destroys acne bacteria and lightens acne scars - A 10-minute massage 2-3 times a week is enough to achieve visible results. Note: the device is not suitable for use on the face.

Red - stimulation

In the second mode you will take care of skin nourishment. The massage using radiofrequency and EMS stimulates blood circulation, thus improving metabolism and increasing cell vitality. Thus stimulated, the body responds with increased collagen production, resulting in more supple and radiant skin.



Use the full power of violet light

Mode 3 comprehensively combines the effects of the first two modes. The combination of two frequencies (red and blue light) provides multi-layered care for best results. Massage is best done after bathing and gently moisturizing the body with lotion. Create a routine that will leave you feeling like a goddess.

Safe use

Everyone's skin is different, so it's a good idea to follow some basic rules that will make using the device a pure pleasure. Observe your sensations and don't massage one batch for more than 10 minutes. After this time, the device will automatically turn off to eliminate the risk of irritation. It is a good idea to start the massage with mode 1 and gradually get used to the operation of the device.

Efficient operation

The device uses a rechargeable battery with a capacity of 1100 mAh, which guarantees up to 40 minutes of operation of the massager. In the package you will find a micro USB cable that will help you replenish energy of the device - it will take you about 3 hours.

Telefone: +351962484153

Email: info@mobilemax.pt

(Chamada para rede móvel nacional)

Included

Massager

USB to micro USB cable

Silicone base

Storage bag

Manual

Manufacturer

Anlan

Model

02-AXTY31-02A

Ultrasonic waveform

3.2 MHz

Light wavelength

Red 620±10nm; Blue 470±10nm; Violet 390 ±10nm.

Vibration motor speed

Approximately 9000 rpm.

Charging port

micro USB

Battery capacity

1100 mAh

Operating time

Approximately 40 minutes (maximum power)

Charging time

Approximately 180 minutes

Input

DC 5V 1A

Voltage

7.4V

Smart timer

10 minutes

Degree of protection

IPX6



Preço:

Antes: € 81.4998

Agora: € 75.01

Telefone: +351962484153

Email: info@mobilemax.pt

(Chamada para rede móvel nacional)

Saúde e beleza, Equipamentos para cuidado corporal