



Vibrating gun massager KiCA Mini C (white) Masażer wibracyjny KiCA Mini C do mięśni (biały)

Vibrating Massager KiCA Mini C for Muscles (White)

Say goodbye to persistent pain and experience your dream relaxation. The KiCA Mini C vibrating massager offers a speed of up to 3000 rpm, effectively reaching deep muscle layers. 3 additional attachments allow for efficient massage of various body parts, and the ability to adjust the intensity level makes it easy to customize the device's performance to specific needs. Another advantage is its weight of only 213 g, allowing you to comfortably perform massages for an extended period without discomfort. Additionally, a single charge provides up to 7 hours of operation.

Ref: 6970078072459

Phone: +351962484153

Email: info@mobilemax.pt

(Call to national mobile network)

A Soothing Massage

The KiCA massager provides a range of motion of up to 7 mm and a pressure force of up to 6Kgf, ensuring intensive and satisfying massages. With a maximum speed of 3000 rpm, it effectively reaches deep tissue layers, offering relief from pain and accelerating muscle recovery. The device also supports a Physical therapy mode in which it alternates between 2 speeds - 2000 rpm and 3000 rpm, making it versatile for various situations. The Mini C model proves reliable for those looking to improve their well-being and seek relief from daily tension.

Convenient LED Indicators

The device is equipped with LED indicators that provide real-time information about the selected massage speed. Blue light indicates



that the Mini C is operating at 2000 rpm (Muscle arousal). Green color indicates 2400 rpm (Relax), and red color indicates 3000 rpm (Deep massage). The last indicator - yellow - signals that the Physical therapy mode has been selected, in which the massager automatically changes speeds, ensuring a circulation massage.

Multifunctional Attachments Included

The package includes a matte metal head and 3 silicone attachments, allowing the KiCA massager to care for every muscle group in the body. This versatile device will satisfy the needs of both professional athletes and individuals in need of relaxation after a long day of work.

Lightweight, Compact Design

The massager weighs less than 213 g, and its size is smaller than a smartphone, allowing you to perform massages for an extended period without straining your wrist. The compact design also enables convenient transport - as a result, the Mini C is the perfect travel companion, providing immediate relief from muscle tension anywhere.

Uninterrupted Massage

Thanks to the efficient 1300 mAh battery, the KiCA Mini C massager can work continuously for up to 7 hours (if the speed does not exceed 2000 rpm and you use the flat head). When you opt for higher massage speeds, you can use the device for up to 4.5 hours before needing to recharge its energy. The included USB-C cable allows you to charge it in about 2 hours, and the LED indicators provide information about the battery status.

Safe to Use

To prevent muscle injuries caused by prolonged massage, the KiCA Mini C massager automatically pauses its operation after more than 10 minutes of continuous massage, emitting an audible signal. If it is not restarted within 5 seconds, it will automatically turn off.

Phone: +351962484153

(Call to national mobile network)

Email: info@mobilemax.pt

Included in the Package:

KiCA Mini C Massager with Flat Head

Silicone Ball Head

Silicone Conical Head

Silicone U-shaped Head

USB-C Charging Cable

Storage Bag

User Manual

∏Manufacturer

∏KiCA

∏Model

∏Mini C

□Range of Motion

∏7 mm

□Pressure Force

□6 Kgf

∏Battery

∏7.4 V, 1300 mAh, 9.62 Wh

□Operating Time

□up to 7 hours (may vary depending on massage intensity and attachment used)

□Charging Time

□about 2 hours

☐Working Voltage Range

∏7.1V-8.4V



□Operating Temperature
□from 0°C to 40°C
□Weight
□about 213 g
□Dimensions
□88.2 x 35 x 128 mm
□Color
□White

Price:

€ 50.00

Phone: +351962484153

(Call to national mobile network)

Email: info@mobilemax.pt

Health & Beauty, Massage Guns