



Facial massager EMS ANLAN 01-ADRY11-001

Facial massager EMS ANLAN 01-ADRY11-001

Ref: 6953156300668

ANLAN EMS Facial Massager

Do you dream of smooth, beautiful skin? With ANLAN's bestselling ultrasonic massager in a new version you will get rid of fine lines, firm your complexion and win the fight against acne! The device offers 5 care modes, 3 levels of temperature and 2 types of light therapy. It also allows you to adjust the intensity of EMS microcurrents in 3 steps and uses vibrations with a frequency of about 10,000 rpm. Its operation is easy with the built-in display.

5 available modes

The ANLAN massager offers as many as 5 modes, which will allow you to comprehensively take care of your skin condition. They include Clean (high temperature, vibrations, positive ions, phototherapy), Cool (low temperature, vibrations, phototherapy), EMS lifting (EMS microcurrents, phototherapy), Moisturizing (high temperature, vibrations, negative ions, phototherapy) and Col Synthesis (low temperature, vibrations, negative ions, phototherapy).

Discover the power of ions

The device offers skin care with positive and negative ions. Negative ions make the skin better absorb nutrients. They make cosmetics more effective. Positive ions enable thorough, effective cleansing. They perfectly remove dirt and makeup residues. As a result the skin becomes clean, moisturized and radiant.

Benefit from light therapy

Find out how many benefits phototherapy can bring. A massager can work on the skin with red or blue light. Red light increases skin elasticity. It firms it and prevents the formation of wrinkles. Blue light eliminates bacteria, disinfects and cleanses. Thus, it can reduce redness and also helps in the treatment of acne.

Smooth, supple skin

Say goodbye to wrinkles and minor imperfections. The device by ANLAN offers EMS microcurrent massage, which speeds up metabolism and makes skin cells renew themselves faster. As a result, the complexion becomes more elastic and smooth. Now you can easily restore your skin's firmness and youthful appearance! The massager also allows you to adjust the intensity of EMS microcurrents in 3 steps.

High and low temperature

The device also offers warm and cool compresses. In either case, it also allows you to adjust the temperature in 3 stages. High temperatures (38°C / 40°C / 42°C) allow the pores to open quickly, allowing nutrients to be better absorbed. Low temperatures (6°C / 8°C / 12°C) tighten pores and allow skincare products to last longer.

Clever design

The ANLAN cordless massager fits comfortably in your hand and only takes about 3 hours to charge. The built-in display provides convenient access to information - for example, you can check the battery charge level. A special stand makes storing the device easy. The massager also has air holes that ensure efficient heat dissipation and increase safety of use.

In the box

Massager x1

Cotton Fix Ring x1

USB charging cable x1

User manual x1

Brand

ANLAN

Model

01-ADRY11-001

Material

ABS + stainless steel

Battery capacity

1000mAh

Temperature range

6-42°C

Vibration frequency

10,000rpm±10%

Charging time

About 3h

Weight

135g

Device dimensions

175x55x67mm

Base size

80x65mm

Modes of operation

Clean, Cool, EMS Lifting, Moisturizing, Col Synthesis

Preço:

€ 45.01

Saúde e beleza, Massage Guns